



**Manufacturer:**  
Chemtronics

**Product Name:**  
Coventry Sealed Polyester Swabs, 0.12" Head - 500 Swabs/Bag

**Manufacturer Part Number:**  
38540

▶ Click here for more details on the Coventry Sealed Polyester Swabs, 0.12" Head - 500 Swabs/Bag

# COVENTRY™ TECHNICAL Data Sheet

Rev. C (8/21) Page 1 of 3

## Sealed Polyester Swabs

**Product# 34860, 38040, 38140, 38540, 36060, 36061I, 36850**

### Product Description

Coventry™ Sealed Fabric Swabs are constructed from either a non-woven or a knit fabric. The head is thermoformed and welded to the handle in a single step using a proprietary process that yields an exceptionally clean swab that will not fall apart. The heads are made from either *polyester knit, a polyester/nylon microfiber knit, or a polyester non-woven fabric*. These swabs are acid and solvent resistant, including acetone.

- Excellent strength
- High absorbency and absorbency rate
- Excellent solvent compatibility
- Does not generate loose fibers or particles
- No adhesives or binders are used in the construction
- Good for tight tolerance work
- Micro-fiber has increased number of microspace to trap more particles, regardless of wet or dry
- Available Aqua-Prime™ processed


### Typical Applications

Coventry™ Sealed Polyester Swabs can be used to:

- Validate biogenetic/pharmaceutical surfaces
- Clean disk drive components
- Clean sensitive surfaces in optics
- Clean aerospace hardware

### Compatibility

**Fabric Heads and Polypropylene Handles** are compatible with most common solvents and dilute or weak acids, as well as ketones such as acetone or methyl ethyl ketone. Aqua Prime™ process available on all polyester swabs.



### Availability

- 34860** 3.2" (8.2 cm) polypropylene handle, 0.11" (2.7 mm) conical microfiber head, 500/bag
- 38040** 2.7" (6.9 cm) polypropylene flexible tip paddle handle, 0.12" (3.0 mm) knit polyester fabric head, 500/ bag
- 38140** 2.7" (6.9 cm) polypropylene flexible tip paddle handle, 0.12" (3.0 mm) knit polyester fabric *thin* 0.07" (1.8 mm) head, 500/ bag
- 38540** 2.7 " (6.9 cm) polypropylene flexible tip paddle handle, 0.12" (3.0 mm) non-woven polyester fabric head, 500/ bag
- 36060** 5.8" (14.7 cm) polypropylene handle, 0.21" (5.6 mm) knit polyester fabric head, 500/bag
- 36061I** 36060 Aqua Prime™ processed, 25/bag
- 36850** 5.8" (14.7 cm) glass filled nylon handle, 0.21" (5.3 mm) knit polyester/nylon microfiber fabric head (swab handle tip slightly tapered), 500/bag

Contact the professionals at Fiber Optic Center for a quote or to get more details.

[focenter.com](http://focenter.com) • 508-992-6464 | (800) 473-4237 • [sales@focenter.com](mailto:sales@focenter.com)

23 Centre Street • New Bedford, MA 02740 USA

Product specifications and data are subject to change without notice. FOC last update 3/5/2026.



**Manufacturer:**  
Chemtronics

**Product Name:**  
Coventry Sealed Polyester Swabs, 0.12" Head - 500 Swabs/Bag

**Manufacturer Part Number:**  
38540








▶ [Click here for more details on the Coventry Sealed Polyester Swabs, 0.12" Head - 500 Swabs/Bag](#)

## TECHNICAL Data Sheet

Rev. C (8/21) Page 2 of 3

### Sealed Polyester Swabs

Product# 34860, 38040, 38140, **38540**, **36060**, 36061, 36850




Contact the professionals at Fiber Optic Center for a quote or to get more details.

[focenter.com](http://focenter.com) • 508-992-6464 | (800) 473-4237 • [sales@focenter.com](mailto:sales@focenter.com)

23 Centre Street • New Bedford, MA 02740 USA

Product specifications and data are subject to change without notice. FOC last update 3/5/2026.



**Manufacturer:**  
Chemtronics

**Product Name:**  
Coventry Sealed Polyester Swabs, 0.12" Head - 500 Swabs/Bag

**Manufacturer Part Number:**  
38540




▶ Click here for more details on the Coventry Sealed Polyester Swabs, 0.12" Head - 500 Swabs/Bag

## TECHNICAL Data Sheet

Rev. C (8/21) Page 3 of 3

### Sealed Polyester Swabs

Product# 34860, 38040, 38140, 38540, 36060, 36061I, 36850



### Sealed Polyester Swabs Test Data

Head Material	Polyester Knit Fabric			
Swabs:	36060	36061I	38040	38140
<b>Particle Counts:</b>				
0.5 - 1 um	685	504	685	685
1.0 - 2 um	189	81	189	189
2.0 - 5 um	58	31	58	58
5.0 - 10 um	30	8	30	30
10 - 25 um	19	9	19	19
Over 25 um	9	2	9	9
TOTAL	990	635	990	990
NVR (ug/swab)	none detected	none detected	none detected	none detected
Chlorides (ug/swab)	none detected	1	none detected	none detected
TOC (ug/swab)	4.2	4.1	4.2	4.2

The swabs were tested by Particle Measurement Technology. All the swabs tests were done using the recommended practices of the Institute of Environmental Science, Swabs and Wipes Working Group.

Test data not available for 38540, 34860 and 36850.

#### Note:

This information is believed to be accurate. It is intended for professional end users having the skills to evaluate and use the data properly. CHEMTRONICS does not guarantee the accuracy of the data and assumes no liability in connection with damages incurred while using it.

Chemtronics® is a registered trademark of Chemtronics. All rights reserved. Coventry™ and Aqua Prime™ are trademarks of Chemtronics. All rights reserved.

Contact the professionals at Fiber Optic Center for a quote or to get more details.

[focenter.com](http://focenter.com) • 508-992-6464 | (800) 473-4237 • [sales@focenter.com](mailto:sales@focenter.com)

23 Centre Street • New Bedford, MA 02740 USA

Product specifications and data are subject to change without notice. FOC last update 3/5/2026.