



**Manufacturer:**  
Branson Ultrasonics Corporation

**Product Name:**  
Perforated Tray for Branson M2800 Ultrasonic Cleaner

**Manufacturer Part Number:**  
100-410-162

▶ [Click here for more details on the Perforated Tray for Branson M2800 Ultrasonic Cleaner](#)

## Brasonic® Ultrasonic Cleaning Bath Accessories



### General Description

Brasonic accessories are designed to maximize the effectiveness of your Brasonic baths. It is important that parts do not rest on the bottom of the bath, as it can inhibit the ultrasonic activity and possibly damage the unit.

A variety of perforated trays, mesh baskets, beakers, and support racks give you the option of direct or indirect cleaning. All five models come with covers to allow the cleaning bath to heat up more quickly.

### Direct Cleaning Method

- Suspends parts using a perforated tray, mesh basket, or a support rack.
- Permits sound waves to penetrate directly to the parts.

### Indirect Cleaning Method

- Suspends parts using a solid tray, or beakers with a beaker cover.
- Ideal for cleaning small parts.
- Discarded soil cleans easily from a beaker or solid tray.
- Indirect method is ideal for sample preparation.

### Additional Features

- Perforated trays, solid insert trays, and wire mesh baskets have built-in handles specially designed to prevent accessories from touching the tank bottom.
- Support racks allow total volume of the tank to be used while keeping parts, fixtures, or containers about 1/2" above the tank bottom.
- Support racks also are useful as countertop drain racks, holding parts above the countertop surface.

**Contact the professionals at Fiber Optic Center for a quote or to get more details.**

[focenter.com](http://focenter.com) • 508-992-6464 | (800) 473-4237 • [sales@focenter.com](mailto:sales@focenter.com)

23 Centre Street • New Bedford, MA 02740 USA

Product specifications and data are subject to change without notice. FOC last update 3/8/2026.



**Manufacturer:**  
Branson Ultrasonics Corporation


**Product Name:**  
Perforated Tray for Branson M2800 Ultrasonic Cleaner

**Manufacturer Part Number:**  
100-410-162

▶ [Click here for more details on the Perforated Tray for Branson M2800 Ultrasonic Cleaner](#)

**Brasonic Accessories** – *Order using part numbers listed under model numbers*

	Ordering Part No.	Model 2800
	<b>Cover</b> (Supplied with Unit)	100-032-517
	<b>Perforated Tray</b>	100-410-162
	<b>Mesh Basket</b>	100-916-334
	<b>Solid Insert Tray</b>	100-410-172
	<b>Beaker Positioning Cover</b>	(2x600 ml) 100-032-518
		(2x250 ml) 100-032-527
	<b>Support Racks</b>	CPN-916-040



**Contact the professionals at Fiber Optic Center for a quote or to get more details.**

[focenter.com](http://focenter.com) • 508-992-6464 | (800) 473-4237 • [sales@focenter.com](mailto:sales@focenter.com)  
23 Centre Street • New Bedford, MA 02740 USA

Product specifications and data are subject to change without notice. FOC last update 3/8/2026.