



### **Manufacturer:**

**Branson Ultrasonics Corporation** 

### **Product Name:**

Beaker Cover (2x250 ml) for Branson M2800 Ultrasonic Cleaner

### **Manufacturer Part Number:**

100-032-527

▶ Click here for more details on the Beaker Cover (2x250 ml) for Branson M2800 Ultrasonic Cleaner

# Bransonic® Ultrasonic Cleaning Bath Accessories

# **General Description**

Bransonic accessories are designed to maximize the effectiveness of your Bransonic baths. It is important that parts do not rest on the bottom of the bath, as it can inhibit the ultrasonic activity and possibly damage the unit.

A variety of perforated trays, mesh baskets, beakers, and support racks give you the option of direct or indirect cleaning. All five models come with covers to allow the cleaning bath to heat up more quickly.

# **Direct Cleaning Method**

- Suspends parts using a perforated tray, mesh basket, or a support rack.
- Permits sound waves to penetrate directly to the parts.

# Indirect Cleaning Method

- Suspends parts using a solid tray, or beakers with a beaker cover.
- Ideal for cleaning small parts.
- Discarded soil cleans easily from a beaker or solid tray.
- Indirect method is ideal for sample preparation.



### **Additional Features**

- Perforated trays, solid insert trays, and wire mesh baskets have built-in handles specially designed to prevent accessories from touching the tank bottom.
- Support racks allow total volume of the tank to be used while keeping parts, fixtures, or containers about ½" above the tank bottom.
- Support racks also are useful as countertop drain racks, holding parts above the countertop surface.





# Manufacturer:

**Branson Ultrasonics Corporation** 

# **Product Name:**

Beaker Cover (2x250 ml) for Branson M2800 Ultrasonic Cleaner

# **Manufacturer Part Number:**

100-032-527

Click here for more details on the Beaker Cover (2x250 ml) for Branson M2800 Ultrasonic Cleaner

# **Bransonic Accessories –** Order using part numbers listed under model numbers

Ordering Part No.	Model 2800
<b>Cover</b> (Supplied with Unit)	100-032-517
Perforated Tray	100-410-162
Mesh Basket	100-916-334
Solid Insert Tray	100-410-172
Beaker Positioning Cover	(2x600 ml) 100-032-518 (2x250 ml) 100-032-527
Support Racks	CPN-916-040
	Perforated Tray  Mesh Basket  Solid Insert Tray  Beaker Positioning Cover  Support

Contact the professionals at Fiber Optic Center for a quote or to get more details.